

Traverse Rescue Stretcher I.V. Pole



Walking beside or behind a stretcher on a narrow steep trail holding an I.V. solution bag in the air can be a challenge. The Traverse Rescue Stretcher durable stainless steel I.V. Pole works with one of three different mounting blocks to attach to the top rail of several Traverses basket stretchers. The I.V. Pole is secured to the mounting block with the use of an attached clevis pin and it is easily removable for storage when not use.



Part No.	Description
15-0117	I.V. Pole (22" and _" stainless steel
15-0118	I.V. Pole Mounting Block for Titan & Titan-Ti Stretchers
15-0119	I.V. Pole Mounting Block for Advantage & Bugaboo Stretchers
15-0120	I.V. Pole Mounting Block for Gazelle Stretchers



Traverse Rescue 755 The Queensway E., Ste. 8, Mississauga, Ont. L4Y 4C5
 Email: info@traverserescue.com Website: www.traverserescue.com
 Phone: 905.279.7866 Toll Free (US & Canada) 1.866.998.7866 Fax: 905.615.1367